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22-year-old female patient.

COMPLAINT

Patient reported hypersensitivity to cold in tooth 16.

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EN CLINICAL/RADIOGRAPHIC EXAM

Presence of caries on the pulp wall of the cavity, under an existing restoration.

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TREATMENT

Removal of the old restoration and carious tissue and re-restoration by associating a flowable microhybrid composite in the cavity bottom and conventional microhybrid composite in the restoration remainder.

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STEP BY STEP | PASO A PASO

1. Initial appearance of the tooth showing poor restoration.

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2. Restoration removal and cavity preparation performed, showing an irregular surface on the pulp wall.

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3. Starting with the adhesive procedures, acid etching was performed (Condac 37) for 15 seconds on dentin and 30 seconds on enamel, followed by abundant washing and gentle drying, allowing the dentin to remain moist.

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
4a and 4b. By using a microapplicator (Cavibrush), the adhesive (Ambar APS) was then applied on the moist surface, by rubbing the first drop of the product for 10 seconds. Then, a new layer of adhesive was applied on the same surface for another 10 seconds. The air jet was then applied for 10 seconds, which contributes to the solvent evaporation and, consequently, to the increase of adhesion. The adhesive was then light cured for 10 seconds.

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5a and 5b. To regulate the pulp wall of the cavity and to reduce the curing tension of the composite, a flowable composite was used (Opallis Flow) at the bottom of the cavity, which was inserted and then light cured for 40 seconds.

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6a to 6f. Oblique increments performed in Opallis DA2 for the innermost layers and EA1 for the most superficial layers.

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// The correct use of adhesives, according to the manufacturer's instructions, and the association of flowable and conventional composites, allow a better distribution of the stress generated during composite light curing, which, in theory, increase the durability of the restorations and decreases the risk of postoperative hypersensitivity. //

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Prof. Dr. Marcos Barceleiro

7a and 7b. Checking the occlusion, finishing and polishing the restoration.

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8. Restoration completed.

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FGM MATERIALS USED | MATERIALES FGM UTILIZADOS

LIGHT-CURING FLOWABLE COMPOSITE
Opallis Flow

LIGHT-CURING COMPOSITE
Opallis

5TH GENERATION COLORLESS
LIGHT-CURING ADHESIVE
Ambar APS

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