CASE HIGHLIGHTS | DESTAQUES DEL CASO

Figure 1: Initial picture showing a large dental destruction involving the vestibular, mesial and occlusal faces of the right inferior second pre-molar. The endodontics treatment was already performed and the tooth was restored in an interim basis.
 Figure 2: Color proofing. Based on the tooth remaining and the adjacent teeth, we defined the colors of the resins that would be used with the aid of Opallis shade guide (FGM). In this image, color proofing for dentin and enamel.
 Figure 3: Tooth remaining with guta-percha exposition of the endodontic filling.
 Figure 4: Choice of the DC2 fiberglass post based on the template of the Whitepost (FGM) system.
 Figure 5: Opening and access with the aid of the Gates Glidden drill.
 Figure 6: Note the x-ray and the angulation of the radicular canal. This step is important avoiding trepanation of the root. Beginning of the preparation of the canal with the Whitepost (FGM) drill.
 Figure 7: Post proof and laboratory pencil marking of the high cut of the post.
 Figure 8: Proof of the cut post in position in the radicular canal, avoiding its exposition to the oral medium.

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EN

Opallis onlay by using universal adhesive – A simple alternative for rehabilitation of posterior teeth

AUTHORS: Carlos Eduardo Francci, Guilherme Saavedra and Alexander Nishida

Female patient, 52.

COMPLAINT

Fractured tooth.

INITIAL APPEARANCE

Extensive crown destruction of tooth 45 involving the vestibular peak, reaching the pulpar chamber.

INTRODUCTION

Cured composite resins in ideal conditions can be used in the office as an alternative for the ceramic onlays and the laboratory composite resin restorations. They present a number of advantages, such as less time needed and the lower cost of preparation since they can be made in the office not requiring special equipment. They are easy and quick to make. On top of that, they have technical advantages such as better polishing, better curing degree, better shade stability and a good marginal sealing. This technique associated with the use of a universal adhesive provides greater convenience and versatility to the professional.

CLINICAL/X-RAY EXAM

Good general oral health, presenting fracture on teeth 45 compromising the pulp and no periapical lesion.

TREATMENT

After the removal of the cavity lesion, access and endodontic treatment, reconstruction with the use of fiberglass post (Whitepost) and the cementation of an indirect restoration made with composite resin (Opallis) was planned. The goal was to preserve as much as possible the healthy dental tissue, using, this way, a minimally invasive and low cost technique.

Figure 9: Ambar Universal (FGM) adhesive system application on dentine and enamel. The selective etching technique was used, only on enamel and crown dentin, being the radicular one not etched.

Figure 10: Solvent evaporation with air jet.

Figure 11: Silanization: A layer of the silano Prosil (FGM) was applied to the post and the post was left to dry naturally for one minute.

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
CASE HIGHLIGHTS | DESTAQUES DEL CASO

Figura 12: One drop of Ambar Universal (FGM) adhesive system was placed on the Cavibrush (FGM) tip for the application on the fiberglass post. (Ambar Universal contains MDP, which provides stable and efficient adhesion on enamel and dentine, as well as on a variety of restorative materials).

Figura 13: Application of the of the cement Allcem (FGM) since the radicular canal base.

Figura 14: Placing the post into position. It must be inserted with pressure to make sure it goes all the way to the end.

Figura 15: Injection of Allcem Core (FGM) around the fiberglass post for the filling of the crown portion and, after that, light curing for 40 seconds.

Figura 16: Dental preparation with Diamond tip FG (#2136) tapered rounded tip.

Figura 17: Hemosthase (FGM), hemostatic gel used to contain gingival bleeding.

Figura 18: Preparation ready for molding.

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Figura 18: Preparo listo para el moldeo.

Figure 19: Composite being accommodated (already over the plaster model) against the preparation with the aid of the rubber-tipped brush.

Figure 20: Aspect of the indirect restoration with the dentine base and interface built.

Figure 21: Finishing of the proximal.

Figure 22: Restoration separated from the model (cementation aerial view).

Figure 23: Application of Ambar Universal (FGM) adhesive system on the piece.

Figure 24: Detail of the application of Ambar Universal (FGM) adhesive system on the tooth.

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
Figure 25: Application of Allcem Core (FGM) on the prepared tooth.

Figure 26: Application of the cement on the piece.

Figure 27: Light curing for 40 seconds.

Figure 28: Fine finishing with Diamond Pro (FGM) sandpaper disks.

Figure 29: Polishing with diamond paste Diamond ACI (FGM) with Diamond Feltre (FGM) in low rotation.

Figure 30: Polishing with Diamond ACII (FGM) paste.

Figure 31: Polishing with final polishing paste Diamond R (FGM) and Diamond Feltre (FGM) in low rotation.

Figure 32: Adequate points of contact after cementation.

Figure 33: Lingual view of the cemented restoration.

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Once again, FGM improves what it was already doing well. Following market trends, FGM launches Ambar Universal, a versatile adhesive with the already known quality of Ambar, with MDP, ethanol solvent, excellent consistency and nanoparticles fillers. Now with the acidic power of a self-etching adhesive, which allows for the selective technique where professionals can etch the enamel with Condac 37 but the dentin is only etched by Ambar Universal itself. Another gift from FGM's 20th anniversary. Congratulations to us, Brazilian dentists, who are very proud of having a Brazilian company that represents us so well in the global scenario.

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