[EN] The importance of predictability to the patient by using the CSD

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26-year-old female patient.

COMPLAINT

Patient was not satisfied with the appearance of her teeth.

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[EN] CLINICAL/RADIOGRAPHIC EXAM

The patient came to the office asking for better aesthetics for her smile. She had just completed an orthodontic treatment that did not please her. There was a visible disharmony of the anterior teeth and morphological alterations of teeth 11 and 21. There was a divergent format of the smile. Healthy periodontium.

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TREATMENT

After the clinical and radiological diagnostics, study models were made to better analyze the case and elaborate a treatment plan. A Custom Smile Design (CSD) plan was made with the purpose of achieving a better predictability of the results. After the mock-up confection and the patient's approval, a whitening treatment was carried out in order to improve the color of the substrate. 10 days after finishing the whitening procedure, aesthetic restorations were carried out with composite resin, bringing back the patient's self-esteem.

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1. Patient's initial picture.

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2. Intra-oral initial image.

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3. Planning according to the CSD achieves high predictability in the aesthetic treatments, since it gives evidence to the aesthetic and functional conflict of the patient and directs the professionals to the occlusal approach to follow.

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4. Waxing according to the measurements obtained in the digital drawing.

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5. Molding of the waxing for making the mock-up.

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6. Mock-up: restoration test to obtain the patient's final corroboration, decide on fine tuning details of the form, size and texture.

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7. Whitening was carried out using the mixed technique: 1 in-office application of Whiteness HP Maxx for 45 minutes and 25 days of at-home whitening with Whiteness Perfect 10%.

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
9. Absolute isolation of the teeth followed by deproteinization of the enamel with NaClO at 5.25%, acid etching with Condac 37 (phosphoric acid at 37%) and Clorhexidina at 2% was applied for 30 seconds, with the purpose of inhibiting the metalloproteinases (MMPs), responsible for adhesive degradation.
10 and 11. Application of the Ambar adhesive.

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12. Light curing for 20 seconds.

13. Through a silicone guide the palatine face of the substrate to be restored is done, in this case, with Opallis EA1.

14. Stratification map: effect composites are used to recreate a visual effect of depth, always associated to the opalescence phenomenon.

15. Finalization: checklist of aesthetic parameters. If necessary, corrections are made. Removal of excess with multi-laminated drills and Diamond Pro disks in decreasing order according to their roughness.

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16. One hint to recreate the superficial micro-texture is to mark with a bicolor marker the teeth and do it with multi-laminated drills to corroborate the similarity of the texturization of each of the teeth to be restored. It is recommended to use abrasive rubber to smooth the texture achieved.
17. For final polishing, Diamond Flex felt disks were used in combination with Diamond Excel polishing paste.

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18 and 19. Final result with hydrated teeth. It is possible to observe the excellent shine of the restorative material as well as its mimetization.

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FGM MATERIALS USED | MATERIALES FGM UTILIZADOS

IN-OFFICE WHITENING GEL
Whiteness HP Maxx

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TAKE-HOME WHITENING GEL
(CARBAMIDE PEROXIDE)
Whiteness Perfect 10%

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LIGHT-CURING NANOHYBRID
COMPOSITE
Opallis

RESINA NANOHÍBRIDA FOTOCURABLE
Opallis

5th GENERATION LIGHT-CURING
ADHESIVE
Ambar

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PLATES FOR TRAYS
**Whiteness Plates for Custom
Dental Trays**

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